

## Activity 11 Pre-Lab — Synthesis of Luminol (Forensic Chemistry)

[Your name] & [partner]  
[MM/DD/YY]

Organic Lab · two-step synthesis + chemiluminescence

■ **Black = write this in your notebook** ■ **Blue = context, don't copy it** Copy the header onto every page. Draw dashed-box structures yourself.

### Purpose

Synthesize **luminol** from 3-nitrophthalic acid in two steps, then test its **chemiluminescence** with a blood-substitute (iron catalyst) + hydrogen peroxide — the glow used in forensics to detect blood.

🔥 **DRAW** — the synthesis (2 steps)

- 1) 3-nitrophthalic acid + hydrazine → **3-nitrophthalhydrazide** + 2 H<sub>2</sub>O in triethylene glycol, heat to ~200 °C (distill off water)
- 2) 3-nitrophthalhydrazide + Na<sub>2</sub>S<sub>2</sub>O<sub>4</sub> (in NaOH) → **luminol** (5-amino-2,3-dihydrophthalazine-1,4-dione) reduces –NO<sub>2</sub> → –NH<sub>2</sub>; then acidify w/ acetic acid to precipitate

**CONTEXT** ▶ Directions also ask you to **draw the mechanism** for the boxed portion of the reaction — leave room for it.

### Chemical Reagents & Safety Table

**CONTEXT** ▶ Required (8): 3-nitrophthalic acid, hydrazine, triethylene glycol, aq. NaOH, sodium hydrosulfite, glacial acetic acid, potassium ferricyanide, 3% H<sub>2</sub>O<sub>2</sub>. Typical SDS values — cross-check the SDS tile. **Hydrazine is toxic/carcinogenic** — handle carefully.

| Chemical                         | Role                                        | MW     | Amount             | Physical properties                 | Hazards (GHS)                                                                                  |
|----------------------------------|---------------------------------------------|--------|--------------------|-------------------------------------|------------------------------------------------------------------------------------------------|
| 3-Nitrophthalic acid             | Starting material                           | 211.13 | ~1 g               | Off-white / pale-yellow solid       | H302 harmful if swallowed · H315 · H319 · H335                                                 |
| Hydrazine (8% aq)                | Forms the hydrazide ring                    | 32.05  | 2 mL               | Colorless aqueous liquid            | H301+H311+H331 toxic · H314 corrosive · H317 sensitizer · H350 may cause cancer · H410 aquatic |
| Triethylene glycol               | High-boiling solvent                        | 150.17 | 3 mL               | Colorless viscous liquid; bp 285 °C | No significant GHS hazards (high-boiling solvent)                                              |
| Sodium hydrosulfite (dithionite) | Reduces –NO <sub>2</sub> → –NH <sub>2</sub> | 174.11 | as directed        | White powder                        | H251 self-heating (may ignite) · H302 · EUH031 toxic gas w/ acid                               |
| Sodium hydroxide (aq)            | Base for reduction step                     | 40.00  | as directed        | Colorless aqueous solution          | H290 · H314 severe skin burns & eye damage                                                     |
| Glacial acetic acid              | Acidify → precipitate luminol               | 60.05  | as directed        | Colorless liquid; bp 118 °C         | H226 flammable · H314 severe skin burns & eye damage                                           |
| Potassium ferricyanide           | Catalyst (blood/iron mimic)                 | 329.24 | Part 2 (glow test) | Red crystalline solid               | Low hazard · EUH032 — do NOT mix with acid (releases toxic gas)                                |
| Hydrogen peroxide (3%)           | Oxidant for the glow                        | 34.01  | Part 2 (glow test) | Colorless liquid                    | H319 eye irritation (dilute)                                                                   |

**GRAB FIRST:** 50-mL round-bottom (no stir bar) · heating mantle + Variac + raised lab jack · simple-distillation head + **high-temp thermometer ( $\geq 200^\circ\text{C}$ )** + condenser + receiver · boiling stones · Erlenmeyer + steam bath (heat 15 mL water) · ice bath · Büchner + filter flask + vacuum · 3-mL syringe (hydrazine) + 25-mL Erlenmeyer.

### Procedure Plan (left = plan, right = fill in during lab)

#### Plan (write before lab)

##### Reaction 1 — nitro-luminol

- Heat **15 mL water** on a steam bath (boiling stones) — needed later (step 5).
- Put **~1 g 3-nitrophthalic acid** in a **50-mL RBF** (no stir bar); clamp; lower into a heating mantle on a raised lab jack.
- At the hydrazine hood (nitrile gloves) draw **2 mL 8% hydrazine** into a syringe; add it to the flask. **Hydrazine is toxic + a carcinogen — gloves, no skin/eye contact, fresh gloves after.**
- Add a few boiling stones; gently heat (Variac 75–80%) until the acid dissolves. Add **3 mL triethylene glycol**.
- Build a **simple distillation** (high-temp thermometer **submerged** in the solution, not touching bottom). Instructor checks → Variac 80%, heat until the vapor/solution reads **200 °C** (water distills off) = reaction done. **Don't adjust thermometer while hot — remove, tweak, replace.**
- Remove from heat, cool to ~100 °C, add the **15 mL hot water**; cool in an ice bath; **vacuum-filter** the light-yellow nitro compound (no need to dry).

##### Reaction 2 — reduce to luminol

- Dissolve the nitro compound in **NaOH**, add **sodium hydrosulfite** (heat as directed) to reduce  $-\text{NO}_2 \rightarrow -\text{NH}_2$ ; then **acidify with acetic acid** to precipitate **luminol**; vacuum-filter, weigh → % yield.

##### Part 2 — chemiluminescence

- Make the luminol glow: basic luminol + **3%  $\text{H}_2\text{O}_2$**  +  **$\text{K}_3\text{Fe}(\text{CN})_6$**  catalyst → blue light. Observe/record. **Keep  $\text{K}_3\text{Fe}(\text{CN})_6$  away from acid.**

🔥 **LEAVE ROOM**

for the required mechanism drawing of the boxed reaction step.

#### Experiment Notes (in lab)

Mass 3-nitrophthalic acid: \_\_\_\_\_ g

Rxn 1 — max temp reached / time: \_\_\_\_\_

Nitro compound appearance:

\_\_\_\_\_

Rxn 2 — color changes on reduction/acidify:


\_\_\_\_\_

Luminol mass / % yield: \_\_\_\_\_

Glow test — color/brightness/duration:

\_\_\_\_\_

### Setup — draw these



simple distillation (Rxn 1)

#### CONTEXT (don't copy) ▶

**Why distill to 200 °C?** the hydrazide ring only forms at high temp; triethylene glycol (bp 285 °C) lets you get there, and boiling off the water pushes the equilibrium to product.

**Why does luminol glow?**  $\text{H}_2\text{O}_2$  (with the iron catalyst) oxidizes luminol to an excited 3-aminophthalate that drops to ground state and emits blue light — the same reaction forensics uses to find blood (iron in hemoglobin is the catalyst).

### Clean-up & Conclusions

- Hydrazine rinses → labeled hydrazine waste; distillate water → aqueous waste; keep ferricyanide from acid; organic/aqueous to correct waste; bench task.
- Conclusions** (~4 sentences **after lab**): overall % yield of luminol; what each step did (ring formation, then  $-\text{NO}_2 \rightarrow -\text{NH}_2$  reduction); what the glow test demonstrated ( $\text{H}_2\text{O}_2$  + iron catalyst → chemiluminescence); a loss/error source.