

## Activity 7 Pre-Lab — Grignard: Synthesis of Triphenylmethanol

[Your name] & [partner]  
[MM/DD/YY]

Organic Lab · Grignard reagent + addition to a ketone


■ Black = write this in your notebook ■ Blue = context, don't copy it Copy the header onto every page. Draw dashed-box structures yourself.

### Purpose


Make the **Grignard reagent phenylmagnesium bromide** from bromobenzene + Mg, then use it as a **carbon nucleophile** that adds to the carbonyl of **benzophenone**. Acidic work-up gives **triphenylmethanol**, purified by recrystallization.

🔥 DRAW — reproduce both steps in your notebook (from the procedure)

#### Part 1. Grignard reagent synthesis



#### Part 2. Carbonyl addition reaction



Part 1 makes PhMgBr (dry Et<sub>2</sub>O, I<sub>2</sub> initiator); Part 2: PhMgBr adds to benzophenone, then dilute H<sub>2</sub>SO<sub>4</sub> work-up protonates the alkoxide → triphenylmethanol.

CONTEXT ▶ ALL glassware must be bone-dry — water quenches (destroys) the Grignard. Don't wash glassware the day of lab.

### Chemical Reagents & Safety Table

CONTEXT ▶ Required: magnesium, iodine, diethyl ether, bromobenzene, benzophenone, sulfuric acid. Typical SDS values — cross-check the SDS tile.

| Chemical                     | Role                                 | MW     | Amount used                                   | Physical properties                                        | GHS pictograms                                                                                                                                                                                                                                                                                         |
|------------------------------|--------------------------------------|--------|-----------------------------------------------|------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Magnesium (turnings)         | Makes the Grignard                   | 24.31  | 0.5 g (~0.021 mol)                            | Silvery metal solid                                        | <br>flammable                                                                                                                                                                                                     |
| Iodine                       | Initiator (activates Mg)             | 253.81 | 1 small crystal                               | Purple-black solid; sublimes                               | <br>irritant <br>health <br>environ.    |
| Diethyl ether (anhydrous)    | Dry solvent                          | 74.12  | 12 + 5 + 30 mL                                | Colorless liquid; bp 34.6 °C; d 0.71 g/mL; flash pt -45 °C | <br>flammable <br>irritant                                                                                                   |
| Bromobenzene                 | Grignard precursor (limiting)        | 157.01 | 2.0 mL (~3.0 g, ~0.019 mol)                   | Colorless liquid; bp 156 °C; d 1.50 g/mL; flash pt 51 °C   | <br>flammable <br>irritant <br>environ. |
| Benzophenone                 | Electrophile (ketone)                | 182.22 | 3.2 g (~0.018 mol) in 30 mL Et <sub>2</sub> O | White solid; mp 48 °C                                      | <br>irritant <br>health <br>environ.    |
| Sulfuric acid (conc. → dil.) | Acidic work-up (protonates alkoxide) | 98.08  | a few mL conc. → into 7 mL water              | Colorless viscous liquid; d 1.84 g/mL                      | <br>corrosive                                                                                                                                                                                                     |

FYI ▶ also used: anhydrous CaCl<sub>2</sub> (drying tube), methanol (recrystallization solvent).

**GRAB FIRST (all DRY):** 250-mL round-bottom · stir bar + stir plate · Claisen adapter · condenser + thermometer adapter · drying tube + cotton + anhydrous  $\text{CaCl}_2$  · 125-mL dropping funnel · iron rings + clamps · Tygon tubing · steam bath (Part 2) · Büchner + filter flask (Part 3).

### Procedure Plan (left = plan, right = fill in during lab)

#### Plan (write before lab)

##### Part 1 — Make $\text{PhMgBr}$

1. **0.5 g Mg turnings** + stir bar in a **dry 250-mL RBF**; stir vigorously 5–10 min (stopper on). *Grinding strips the  $\text{MgO}$  coat.*
2. Pack a **drying tube**: cotton plug → anhydrous  $\text{CaCl}_2$  (not tight) → cotton.
3. Instructor drops in **1 small  $\text{I}_2$  crystal** (initiator).
4. Build the **reflux apparatus**: Claisen on the RBF → straight arm gets condenser + thermometer adapter + drying tube; curved arm gets the **125-mL dropping funnel**. Turn condenser water to a steady drip.
5. Add solutions via the funnel **off the apparatus** each time (stopcock closed while filling): first **12 mL anhydrous  $\text{Et}_2\text{O}$**  → run onto the Mg. Note color changes.
6. Next add **2.0 mL bromobenzene + 5 mL  $\text{Et}_2\text{O}$**  dropwise; the mix should self-reflux (cloudy/grey =  $\text{PhMgBr}$  forming). *If it won't start: warm gently / crush Mg with a stir rod.*

##### Part 2 — Add benzophenone

7. Add **3.2 g benzophenone in 30 mL  $\text{Et}_2\text{O}$**  dropwise via the funnel; reflux (steam bath) as directed.

##### Part 3 — Work-up + purify

8. Prepare cold dilute acid: add conc.  $\text{H}_2\text{SO}_4$  to **7 mL DI water** (acid to water!). Remove drying tube; add all the cold acid to quench the alkoxide →  $\text{Ph}_3\text{C-OH}$ . *Conc.  $\text{H}_2\text{SO}_4$  is extremely corrosive.*
9. Isolate crude solid; **recrystallize from methanol** (see Activity 3). Weigh → **% yield**; take mp / IR.

#### ✍️ WRITE

limiting reagent = bromobenzene (0.019 mol) → theoretical  $\text{Ph}_3\text{C-OH}$ ; % yield = actual ÷ theoretical × 100

#### Experiment Notes (in lab)


| measurement                         | value |
|-------------------------------------|-------|
| Mass Mg                             | g     |
| Vol bromobenzene                    | mL    |
| Crude $\text{Ph}_3\text{C-OH}$ mass | g     |
| Recrystallized mass                 | g     |
| mp of product                       | °C    |
| % yield                             |       |

#### Colors on adding $\text{Et}_2\text{O}$ / bromobenzene:

#### Signs Grignard formed (reflux? cloudy grey?):

#### On adding benzophenone / work-up:

### Setup — draw this



reflux apparatus (drying tube · condenser · dropping funnel · Claisen)

#### CONTEXT (don't copy) ▶

**Why dry / drying tube?**  $\text{PhMgBr}$  is a strong base — any water protonates it to benzene and it's dead. The  $\text{CaCl}_2$  tube keeps humid air out.

**Why  $\text{I}_2$ ?** a tiny bit of iodine cleans the Mg surface so the reaction starts.

**Why add acid last?** the first product is an alkoxide salt;  $\text{H}_3\text{O}^+$  protonates it to the neutral alcohol (triphenylmethanol).

### Clean-up & Conclusions

- Ether & organic filtrates → organic waste; acid aqueous → neutralize/aqueous waste; return glassware; bench task.
- **Conclusions** (~4 sentences [after lab](#)): % yield + why (moisture killed some Grignard, transfer losses); what mp/IR confirm about purity/identity; why anhydrous conditions mattered.